

Now
FDA Cleared

ePlex[®]

The ePlex[®] Respiratory Pathogen (RP) Panel

detects and identifies the most common viral and bacterial organisms associated with upper respiratory infection.

Rapid Results for Improved Patient Care



<2 minutes of hands-on time

1

Load sample



2

Load cartridge



3

Report results



To learn more about the ePlex RP Panel, contact us at
info@genmarkdx.com



GenMark Diagnostics, Inc.
5964 La Place Court
Carlsbad, CA 92008
USA

+1 760 448 4300
info@genmarkdx.com
genmarkdx.com

GenMark Diagnostics Europe GmbH
Chollerstrasse 4
6300 Zug
Switzerland

+41 41 530 65 26
info@genmarkdx.com
genmarkdx.com

Viral Targets

Adenovirus
Coronavirus (229E, HKU1, NL63, OC43)
Human Metapneumovirus
Human Rhinovirus/Enterovirus
Influenza A
Influenza A H1
Influenza A H1-2009
Influenza A H3
Influenza B
Parainfluenza 1
Parainfluenza 2
Parainfluenza 3
Parainfluenza 4
Respiratory Syncytial Virus A
Respiratory Syncytial Virus B

Bacterial Targets

Chlamydia pneumoniae
Mycoplasma pneumoniae

